 **BAUERFEIND AG**
Triebeser Straße 16
07937 Zeulenroda-Triebes
Germany
P +49 (0) 36628 66-40 00
F +49 (0) 36628 66-44 99
E info@bauerfeind.com
BAUERFEIND.COM


ErgoPad® weightflex 2

Gebrauchsanweisung
Instruction for use
Notice d'utilisation
Instrucciones de uso


**MADE
IN
GERMANY**


a

b

c

d

(de) deutsch.....3	(hr) hrvatski.....28
(en) english.....5	(sr) srpski.....30
(fr) français.....7	(sl) slovenski.....32
(nl) nederlands.....9	(ro) romanian.....34
(it) italiano.....11	(tr) türkçe.....36
(es) español.....13	(ru) русский.....38
(sv) svenska.....15	(et) eesti.....40
(no) norsk.....16	(lv) latviešu.....42
(fi) suomi.....18	(lt) lietuvių.....44
(da) dansk.....20	(uk) українська.....46
(pl) polski.....22	(ja) 日本語.....48
(cs) český.....24	(ko) 한국어.....50
(sk) slovensky.....26	(zh) 简体中文.....52

(de) deutsch

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en english

These instructions for use can also be downloaded from:

www.bauerfeind.com/downloads

Dear Customer,

Please read through the instructions for use carefully and direct any questions to your physician, medical retailer, or our customer service department at telephone number 0800 0 01 05 30 (calls from Germany are free of charge) or e-mail info@bauerfeind.com.

Intended use

Orthopedic orthoses are functional orthoses for the feet and are suitable for use in treating and preventing foot complaints.

Indications

ErgoPad weightflex 2 foot orthoses provide effective support for the treatment of foot deformities and complaints, such as arch decreased / splay foot. ErgoPad weightflex 2 foot orthoses also encourage the natural gait pattern and therefore increase the wearing comfort of your everyday and business shoes.

Risks of using this product

⚠ Caution: Before using / applying a foot orthosis, you should consult a physician or specialist* on the risks involved. Only a trained specialist may carry out the fitting of the orthosis and provide instruction.

If the available internal volume of the shoe is not sufficient for the foot and the selected foot orthosis, or if the foot orthoses are worn in shoes of an incorrect length and width (or

shoes with inadequate lateral guidance in the case of open shoe types), problems relating to poor fit may occur (e.g. a feeling of constriction or pressure, rubbing).

If the feet are positioned on the edges of the foot orthosis, blisters and calluses may form on the soles of the feet. This may also occur if top covers or cushioning material develop holes, creases, or break due to mechanical overloading. You must stop wearing foot orthoses if they are damaged in this way!

Contraindications

Diabetics and people with rheumatoid foot deformities must seek medical advice before using this product to ensure optimum treatment.

Application instructions

Shoes suitable for wearing with foot orthoses must have sufficient internal volume for both the foot and the selected foot orthosis (adequate shoe length / width and, if necessary, a removable, sufficiently thick cushioned sole, as well as lateral guidance for the foot orthoses in the case of open shoe types). An uncomfortable fit that constricts the foot must always be avoided. **a**

When selecting the correct size of foot orthosis, the packaging information serves merely as a guide. Just like correctly fitting shoes, foot orthoses must be at least 1.0 cm longer than the foot, when the foot is in a weight-bearing, standing position. Please check the length allowance for both feet, as they may differ in length. In case of doubt, base your size selection on the longer foot. **b**

The foot orthosis must be positioned straight in the shoe and have as much surface contact as possible with the shoe insole. Lateral contact between the foot orthosis and the shoe upper helps reduce the risk of pressure marks or blistering on the skin of the soles of the feet. **c**

Foot orthoses must always be worn as a pair! Any exceptions to this must be stipulated by the treating physician.

Please regularly check your feet for any pressure marks, redness or other changes caused by the foot orthosis. If you have any questions, please contact a specialist or your physician.

Adjustment phase: As with many medical products, the ErgoPad weightflex 2 may feel somewhat strange to wear at first. We therefore recommend that you start by wearing it for 1–2 hours at a time and slowly increase the wearing period each time until the foot orthosis feels comfortable.

The running-in period should not last much longer than a week. If you still find it uncomfortable to wear after this period, please contact a specialist.

Foot orthoses must be checked by a specialist after they have been worn for approximately one year to ensure that they are still functioning properly.

Cleaning instructions

The ErgoPad weightflex 2 can be cleaned by hand using soap and lukewarm water. Then leave the foot orthoses to air dry. Do not use any external heat sources (e.g. heaters, hairdryers or direct sunlight). The foot orthoses cannot be washed in a washing machine.

Mild disinfectants suitable for human skin may be applied sparingly (sprayed on). However, this may lead to discoloration of the top cover materials.

Part of the body this product is used for

The foot orthoses must be worn in suitable footwear and with socks or stockings on your feet.

Notes on reuse

The ErgoPad weightflex 2 foot orthoses are designed to be worn repeatedly by one user only.

Warranty

The legal provisions of the country in which the product was purchased apply. Country-specific warranty regulations that govern the relationship between the retailer and the buyer also apply where relevant. If you believe that a claim may be made under the warranty, please first contact the retailer from which you purchased the product directly.

Please do not make any changes to the product yourself. This particularly applies to individual adjustments made by specialists. Please follow our instructions for use and care. These are based on our experience gained over many years and ensure the functionality of our medically effective products for a long time. Only aids that function optimally can help you. Failure to follow the instructions may also restrict the warranty.

Discoloration, in particular fading of the top covers of the foot orthoses, may occur during their use. This is dependent on a number of factors, such as moisture levels and sock / stocking materials. This is not covered by warranty.

Disposal

You can dispose of the product in accordance with national legal requirements.

SPECIALISTS

Modification instructions

Fitting:

For an optimum fit, there should be a length allowance of approximately 1 cm between the foot and the foot orthosis and / or shoe in a standing position. **b**

If the foot orthoses are not the right size for your shoes, you can cut them down or reshape them using thermoforming. Please follow the instructions below.

Cutting:

1. Take the inlays out of your shoes and place them on your ErgoPad weightflex 2.
2. Align the heel sections so that they are positioned precisely on top of one another and draw around the inlay on your ErgoPad weightflex 2.
3. Cut along the outline you have drawn. (Make sure you do not cut the foot orthosis too small. We recommend cutting off less material than you need to at first and trying on the foot orthosis before carrying on cutting).
4. Place the ErgoPad weightflex 2 foot orthoses in your shoes (orthoses for left or right feet are marked as such on the back).

Thermoforming **d**

In exceptional cases and having sought specialist advice, it may be necessary to use thermoforming to adapt the foot orthoses to suit the shape of your feet.

1. After cutting down your ErgoPad weightflex 2, preheat your oven to 90 °C / 200 °F.
2. Put the foot orthoses (with the top covers facing downward on a clean surface) in the oven for approximately 2 minutes.
⚠ Risk of burning! Wear protective gloves when taking the foot orthoses out of the oven.
3. As soon as you have taken the foot orthoses out of the oven, place them immediately in your shoes.
4. Put on your shoes straight away (make sure you are wearing socks or stocking).
5. Then stand for 2 minutes with your feet shoulder-width apart and pointing forward.

Your foot orthoses should now be adapted to the shape of your feet.

Version: 2018-04

* A specialist is any person who is authorized according to the state regulations that are relevant to you for fitting and instruction in the use of orthopedic orthoses.

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Bauerfeind Middle East FZ-LLC
Bauerfeind USA, Inc.