

ORDER FORM

VenoTrain® angioflow

BAUERFEIND AG
 T +49(0)800 001 05 10
 F +49(0)800 001 05 15

Delivery address (if different)

Company

Street, building number

Zip code, city

Date, signature


--	--	--	--	--	--	--	--

Customer number

Patient name

Recommendation: The best way to measure and order the VenoTrain angioflow is with the Bodytronic. This order form serves only as a workaround. All measuring points are required for a perfect fit.

Design: AD, Ccl 1, gFs, caramel color, with or without inspection opening

left	Circumferences	right	Lengths
cD		cD	ℓD
cC0		cC0	ℓC0
cC		cC	ℓC
cB2		cB2	ℓB2
cB1		cB1	ℓB1
cB0		cB0	ℓB0
cB		cB	ℓB
cY	cY		
cA	cA		
Foot length with toe			

- Single Pair
 with inspection opening

All sizes in cm. The purchaser is responsible for ensuring that all size information is correct.

- Measuring point A** – Forefoot along the metatarsophalangeal joint
- Measuring point Y** – Over heel and instep in dorsiflexion with tensed tendons
- Measuring point B** – Narrowest point of the area above the ankle
- Measuring point B0** – Halfway point between measuring points B and B1

- Measuring point B1** – Height of top of calf
- Measuring point B2** – Halfway point between measuring points B1 and C
- Measuring point C** – Largest calf circumference
- Measuring point C0** – Halfway point between measuring points C and D
- Measuring point D** – Narrowest point below the hollow of the knee
- Length of foot** – Length from heel to longest toe